




CD40 Antibody

Product Code	CSB-PA004936LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25942
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 5 protein (21-193AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 5 (B-cell surface antigen CD40) (Bp50) (CD40L receptor) (CDw40) (CD antigen CD40), CD40, TNFRSF5
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD40

Image



Western Blot

Positive WB detected in: Mouse lung tissue,

Mouse stomach tissue


All lanes: TNFRSF5 antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31, 23 kDa

Observed band size: 31 kDa



Immunofluorescent analysis of HeLa cells using CSB-PA004936LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.