



CD40 Antibody, Biotin conjugated

Product Code	CSB-PA004936LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25942
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 5 protein (21-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 5 (B-cell surface antigen CD40) (Bp50) (CD40L receptor) (CDw40) (CD antigen CD40), CD40, TNFRSF5
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD40
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.