




# CDC25B Antibody

<b>Product Code</b>	CSB-PA004995LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P30305
<b>Immunogen</b>	Recombinant Human M-phase inducer phosphatase 2 protein (1-195AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC, IF; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200, IF:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	M-phase inducer phosphatase 2 (EC 3.1.3.48) (Dual specificity phosphatase Cdc25B), CDC25B, CDC25HU2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	CDC25B

## Image



### Western Blot

Positive WB detected in: Jurkat whole cell lysate,

Raji whole cell lysate


All lanes: CDC25B antibody at 3µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 65, 64, 61, 68 kDa

Observed band size: 65 kDa



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA004995LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA004995LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.