



CDKN2A Antibody, Biotin conjugated

Product Code	CSB-PA005089LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42771
Immunogen	Recombinant Human Cyclin-dependent kinase inhibitor 2A protein (1-86AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cyclin-dependent kinase inhibitor 2A (Cyclin-dependent kinase 4 inhibitor A) (CDK4I) (Multiple tumor suppressor 1) (MTS-1) (p16-INK4a) (p16-INK4) (p16INK4A), CDKN2A, CDKN2 MTS1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CDKN2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.