



CHAT Antibody, HRP conjugated

Product Code	CSB-PA005314LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28329
Immunogen	Recombinant Human Choline O-acetyltransferase protein (354-592AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Choline O-acetyltransferase (CHOACTase) (ChAT) (Choline acetylase) (EC 2.3.1.6), CHAT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CHAT