



CX3CR1 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA006236LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P49238 |
| Immunogen | Recombinant Human CX3C chemokine receptor 1 protein (298-351AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | CX3C chemokine receptor 1 (C-X3-C CKR-1) (CX3CR1) (Beta chemokine receptor-like 1) (CMK-BRL-1) (CMK-BRL1) (Fractalkine receptor) (G-protein coupled receptor 13) (V28), CX3CR1, CMKBRL1 GPR13 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Microbiology |
| Target Names | CX3CR1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |