



CX3CR1 Antibody, Biotin conjugated

Product Code	CSB-PA006236LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49238
Immunogen	Recombinant Human CX3C chemokine receptor 1 protein (298-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CX3C chemokine receptor 1 (C-X3-C CKR-1) (CX3CR1) (Beta chemokine receptor-like 1) (CMK-BRL-1) (CMK-BRL1) (Fractalkine receptor) (G-protein coupled receptor 13) (V28), CX3CR1, CMKBRL1 GPR13
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	CX3CR1