



CX3CR1 Antibody, Biotin conjugated

Product Code	CSB-PA006236LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49238
Immunogen	Recombinant Human CX3C chemokine receptor 1 protein (298-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CX3C chemokine receptor 1 (C-X3-C CKR-1) (CX3CR1) (Beta chemokine receptor-like 1) (CMK-BRL-1) (CMK-BRL1) (Fractalkine receptor) (G-protein coupled receptor 13) (V28), CX3CR1, CMKBRL1 GPR13
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	CX3CR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.