



CYP11B1 Antibody, Biotin conjugated

Product Code	CSB-PA006390LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15538
Immunogen	Recombinant Human Cytochrome P450 11B1, mitochondrial protein (88-323AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 11B1, mitochondrial (CYPXIB1) (Cytochrome P-450c11) (Cytochrome P450C11) (Steroid 11-beta-hydroxylase) (EC 1.14.15.4), CYP11B1, S11BH
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP11B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.