



EFEMP1 Antibody, Biotin conjugated

Product Code	CSB-PA007450LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12805
Immunogen	Recombinant Human EGF-containing fibulin-like extracellular matrix protein 1 protein (27-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EGF-containing fibulin-like extracellular matrix protein 1 (Extracellular protein S1-5) (Fibrillin-like protein) (Fibulin-3) (FIBL-3), EFEMP1, FBLN3 FBNL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EFEMP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.