



EFEMP2 Antibody, FITC conjugated

Product Code	CSB-PA007451LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95967
Immunogen	Recombinant Human EGF-containing fibulin-like extracellular matrix protein 2 protein (26-165AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EGF-containing fibulin-like extracellular matrix protein 2 (Fibulin-4) (FIBL-4) (Protein UPH1), EFEMP2, FBLN4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EFEMP2