



EPHA1 Antibody, HRP conjugated

Product Code	CSB-PA007720LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21709
Immunogen	Recombinant Human Ephrin type-A receptor 1 protein (296-414AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin type-A receptor 1 (hEpha1) (EC 2.7.10.1) (EPH tyrosine kinase) (EPH tyrosine kinase 1) (Erythropoietin-producing hepatoma receptor) (Tyrosine-protein kinase receptor EPH), EPHA1, EPH EPHT EPHT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EPHA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.