



EPHB6 Antibody, HRP conjugated

Product Code	CSB-PA007733LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15197
Immunogen	Recombinant Human Ephrin type-B receptor 6 protein (410-553AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin type-B receptor 6 (HEP) (Tyrosine-protein kinase-defective receptor EPH-6), EPHB6
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	EPHB6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.