





ERBB4 Antibody

Product Code	CSB-PA007766LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15303
Immunogen	Recombinant Human Receptor tyrosine-protein kinase erbB-4 protein (1027-1204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor tyrosine-protein kinase erbB-4 (EC 2.7.10.1) (Proto-oncogene-like protein c-ErbB-4) (Tyrosine kinase-type cell surface receptor HER4) (p180erbB4) [Cleaved into: ERBB4 intracellular domain (4ICD) (E4ICD) (s80HER4)], ERBB4, HER4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ERBB4

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA007766LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA007766LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.