



FIGN Antibody, FITC conjugated

Product Code	CSB-PA008676LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5HY92
Immunogen	Recombinant Human Fidgetin protein (132-371AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fidgetin, FIGN
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FIGN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.