



FIS1 Antibody, Biotin conjugated

Product Code	CSB-PA008684LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3D6
Immunogen	Recombinant Human Mitochondrial fission 1 protein (1-122AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial fission 1 protein (FIS1 homolog) (hFis1) (Tetratricopeptide repeat protein 11) (TPR repeat protein 11), FIS1, TTC11
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FIS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.