



GJB1 Antibody, Biotin conjugated

Product Code	CSB-PA009451LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08034
Immunogen	Recombinant Human Gap junction beta-1 protein (215-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gap junction beta-1 protein (Connexin-32) (Cx32) (GAP junction 28 kDa liver protein), GJB1, CX32
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GJB1