



GRIA4 Antibody, HRP conjugated

Product Code	CSB-PA009901LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48058
Immunogen	Recombinant Human Glutamate receptor 4 protein (733-883AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor 4 (GluR-4) (GluR4) (AMPA-selective glutamate receptor 4) (GluR-D) (Glutamate receptor ionotropic, AMPA 4) (GluA4), GRIA4, GLUR4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.