



GRIA4 Antibody, Biotin conjugated

Product Code	CSB-PA009901LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48058
Immunogen	Recombinant Human Glutamate receptor 4 protein (733-883AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor 4 (GluR-4) (GluR4) (AMPA-selective glutamate receptor 4) (GluR-D) (Glutamate receptor ionotropic, AMPA 4) (GluA4), GRIA4, GLUR4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIA4