



GSTA1 Antibody, Biotin conjugated

Product Code	CSB-PA009970HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08263
Immunogen	Recombinant Human Glutathione S-transferase A1 protein (2-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase A1 (EC 2.5.1.18) (GST HA subunit 1) (GST class-alpha member 1) (GST-epsilon) (GSTA1-1) (GTH1) [Cleaved into: Glutathione S-transferase A1, N-terminally processed], GSTA1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GSTA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.