



GSTM1 Antibody, FITC conjugated

Product Code	CSB-PA009979LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09488
Immunogen	Recombinant Human Glutathione S-transferase Mu 1 protein (24-156AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase Mu 1 (EC 2.5.1.18) (GST HB subunit 4) (GST class-mu 1) (GSTM1-1) (GSTM1a-1a) (GSTM1b-1b) (GTH4), GSTM1, GST1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GSTM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.