



GUCY2C Antibody, Biotin conjugated

Product Code	CSB-PA010053LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25092
Immunogen	Recombinant Human Heat-stable enterotoxin receptor protein (636-891AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat-stable enterotoxin receptor (STA receptor) (hSTAR) (EC 4.6.1.2) (Guanylyl cyclase C) (GC-C) (Intestinal guanylate cyclase), GUCY2C, GUC2C STAR
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	GUCY2C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.