



HNRNPLL Antibody, FITC conjugated

Product Code	CSB-PA010630LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVV9
Immunogen	Recombinant Human Heterogeneous nuclear ribonucleoprotein L-like protein (1-72AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heterogeneous nuclear ribonucleoprotein L-like (hnRNPLL) (Stromal RNA-regulating factor), HNRNPLL, HNRPLL SRRF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HNRNPLL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.