



HRAS Antibody, HRP conjugated

Product Code	CSB-PA010726LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01112
Immunogen	Recombinant Human GTPase HRas protein (12-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTPase HRas (H-Ras-1) (Ha-Ras) (Transforming protein p21) (c-H-ras) (p21ras) [Cleaved into: GTPase HRas, N-terminally processed], HRAS, HRAS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HRAS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.