



PITX3 Antibody, Biotin conjugated

Product Code	CSB-PA010844LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75364
Immunogen	Recombinant Human Pituitary homeobox 3 protein (10-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pituitary homeobox 3 (Homeobox protein PITX3) (Paired-like homeodomain transcription factor 3), PITX3, PTX3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PITX3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.