



IDE Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA010988LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P14735 |
| Immunogen | Recombinant Human Insulin-degrading enzyme protein (501-800AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Insulin-degrading enzyme (EC 3.4.24.56) (Abeta-degrading protease) (Insulin protease) (Insulinase) (Insulysin), IDE |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | IDE |