



IDE Antibody, Biotin conjugated

Product Code	CSB-PA010988LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14735
Immunogen	Recombinant Human Insulin-degrading enzyme protein (501-800AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin-degrading enzyme (EC 3.4.24.56) (Abeta-degrading protease) (Insulin protease) (Insulinase) (Insulysin), IDE
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	IDE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.