



IGHMBP2 Antibody, Biotin conjugated

Product Code	CSB-PA011172LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38935
Immunogen	Recombinant Human DNA-binding protein SMUBP-2 protein (646-884AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-binding protein SMUBP-2 (EC 3.6.4.12) (EC 3.6.4.13) (ATP-dependent helicase IGHMBP2) (Glial factor 1) (GF-1) (Immunoglobulin mu-binding protein 2), IGHMBP2, SMBP2 SMUBP2
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	IGHMBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.