





IL33 Antibody

Product Code	CSB-PA011656LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95760
Immunogen	Recombinant Human Interleukin-33 protein (123-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-33 (IL-33) (Interleukin-1 family member 11) (IL-1F11) (Nuclear factor from high endothelial venules) (NF-HEV) [Cleaved into: Interleukin-33 (95-270); Interleukin-33 (99-270); Interleukin-33 (109-270)], IL33, C9orf26 IL1F11 NFHEV
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	IL33

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA011656LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA011656LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.