



# IL33 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA011656LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O95760
<b>Immunogen</b>	Recombinant Human Interleukin-33 protein (123-276AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Interleukin-33 (IL-33) (Interleukin-1 family member 11) (IL-1F11) (Nuclear factor from high endothelial venules) (NF-HEV) [Cleaved into: Interleukin-33 (95-270); Interleukin-33 (99-270); Interleukin-33 (109-270)], IL33, C9orf26 IL1F11 NFHEV
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	IL33
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.