



IL33 Antibody, Biotin conjugated

Product Code	CSB-PA011656LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95760
Immunogen	Recombinant Human Interleukin-33 protein (123-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-33 (IL-33) (Interleukin-1 family member 11) (IL-1F11) (Nuclear factor from high endothelial venules) (NF-HEV) [Cleaved into: Interleukin-33 (95-270); Interleukin-33 (99-270); Interleukin-33 (109-270)], IL33, C9orf26 IL1F11 NFHEV
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	IL33