



ITGAM Antibody, Biotin conjugated

Product Code	CSB-PA011876LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11215
Immunogen	Recombinant Human Integrin alpha-M protein (46-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-M (CD11 antigen-like family member B) (CR-3 alpha chain) (Cell surface glycoprotein MAC-1 subunit alpha) (Leukocyte adhesion receptor MO1) (Neutrophil adherence receptor) (CD antigen CD11b), ITGAM, CD11B CR3A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ITGAM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.