



TTI1 Antibody, Biotin conjugated

Product Code	CSB-PA012188LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43156
Immunogen	Recombinant Human TELO2-interacting protein 1 homolog protein (692-858AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TELO2-interacting protein 1 homolog (Protein SMG10), TTI1, KIAA0406 SMG10
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TTI1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.