



KPNA2 Antibody, Biotin conjugated

Product Code	CSB-PA012484LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52292
Immunogen	Recombinant Human Importin subunit alpha-1 protein (2-285AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Importin subunit alpha-1 (Karyopherin subunit alpha-2) (RAG cohort protein 1) (SRP1-alpha), KPNA2, RCH1 SRP1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KPNA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.