





LCE1B Antibody

Product Code	CSB-PA012785LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T7P3
Immunogen	Recombinant Human Late cornified envelope protein 1B protein (1-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Late cornified envelope protein 1B (Late envelope protein 2) (Small proline-rich-like epidermal differentiation complex protein 2A), LCE1B, LEP2 SPRL2A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LCE1B

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA012785LA01HU at dilution of 1:100



Immunofluorescent analysis of A375 cells using CSB-PA012785LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.