



LCE1B Antibody, Biotin conjugated

Product Code	CSB-PA012785LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T7P3
Immunogen	Recombinant Human Late cornified envelope protein 1B protein (1-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Late cornified envelope protein 1B (Late envelope protein 2) (Small proline-rich-like epidermal differentiation complex protein 2A), LCE1B, LEP2 SPRL2A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LCE1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.