



PLIN3 Antibody, FITC conjugated

Product Code	CSB-PA013294LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60664
Immunogen	Recombinant Human Perilipin-3 protein (192-324AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Perilipin-3 (47 kDa mannose 6-phosphate receptor-binding protein) (47 kDa MPR-binding protein) (Cargo selection protein TIP47) (Mannose-6-phosphate receptor-binding protein 1) (Placental protein 17) (PP17), PLIN3, M6PRBP1 TIP47
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLIN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.