




MAK Antibody

Product Code	CSB-PA013369LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20794
Immunogen	Recombinant Human Serine/threonine-protein kinase MAK protein (200-440AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase MAK (EC 2.7.11.1) (Male germ cell-associated kinase), MAK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAK

Image



Western Blot

Positive WB detected in: Mouse brain tissue


All lanes: MAK antibody at 2.8µg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 71, 74, 67 kDa

Observed band size: 71 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA013369LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA013369LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)