



MAK Antibody, FITC conjugated

Product Code	CSB-PA013369LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20794
Immunogen	Recombinant Human Serine/threonine-protein kinase MAK protein (200-440AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase MAK (EC 2.7.11.1) (Male germ cell-associated kinase), MAK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.