



MAP1LC3B2 Antibody, Biotin conjugated

Product Code	CSB-PA013404LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NCE7
Immunogen	Recombinant Human Microtubule-associated proteins 1A/1B light chain 3 beta 2 protein (1-47AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated proteins 1A/1B light chain 3 beta 2 (Microtubule-associated proteins 1A/1B light chain 3B-like), MAP1LC3B2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAP1LC3B2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.