





MAP3K5 Antibody

Product Code	CSB-PA013427LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99683
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 5 protein (363-431AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 5 (EC 2.7.11.25) (Apoptosis signal-regulating kinase 1) (ASK-1) (MAPK/ERK kinase kinase 5) (MEK kinase 5) (MEKK 5), MAP3K5, ASK1 MAPKKK5 MEKK5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAP3K5

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA013427LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA013427LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.