



MAP3K5 Antibody, HRP conjugated

Product Code	CSB-PA013427LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99683
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 5 protein (363-431AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase kinase kinase 5 (EC 2.7.11.25) (Apoptosis signal-regulating kinase 1) (ASK-1) (MAPK/ERK kinase kinase 5) (MEK kinase 5) (MEKK 5), MAP3K5, ASK1 MAPKKK5 MEKK5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MAP3K5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.