



MAPK10 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA013449LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P53779 |
| Immunogen | Recombinant Human Mitogen-activated protein kinase 10 protein (133-316AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Mitogen-activated protein kinase 10 (MAP kinase 10) (MAPK 10) (EC 2.7.11.24) (MAP kinase p49 3F12) (Stress-activated protein kinase 1b) (SAPK1b) (Stress-activated protein kinase JNK3) (c-Jun N-terminal kinase 3), MAPK10, JNK3 JNK3A PRKM10 SAPK1B |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | MAPK10 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |