



YWHAZ Antibody, HRP conjugated

Product Code	CSB-PA01459B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P63104
Immunogen	Recombinant Human 14-3-3 protein zeta/delta protein (133-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	14-3-3 protein zeta/delta (Protein kinase C inhibitor protein 1) (KCIP-1), YWHAZ
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	YWHAZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.