



MLYCD Antibody, FITC conjugated

Product Code	CSB-PA014646LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95822
Immunogen	Recombinant Human Malonyl-CoA decarboxylase, mitochondrial protein (40-493AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Malonyl-CoA decarboxylase, mitochondrial (MCD) (EC 4.1.1.9), MLYCD
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MLYCD