



MLYCD Antibody, Biotin conjugated

Product Code	CSB-PA014646LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95822
Immunogen	Recombinant Human Malonyl-CoA decarboxylase, mitochondrial protein (40-493AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Malonyl-CoA decarboxylase, mitochondrial (MCD) (EC 4.1.1.9), MLYCD
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MLYCD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.