



MPZL1 Antibody, Biotin conjugated

Product Code	CSB-PA014775LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95297
Immunogen	Recombinant Human Myelin protein zero-like protein 1 protein (195-258AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin protein zero-like protein 1 (Protein zero-related), MPZL1, PZR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MPZL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.