



MSR1 Antibody, Biotin conjugated

Product Code	CSB-PA015050LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21757
Immunogen	Recombinant Human Macrophage scavenger receptor types I and II protein (271-446AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrophage scavenger receptor types I and II (Macrophage acetylated LDL receptor I and II) (Scavenger receptor class A member 1) (CD antigen CD204), MSR1, SCARA1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MSR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.