



MSR1 Antibody, Biotin conjugated

Product Code	CSB-PA015050LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21757
Immunogen	Recombinant Human Macrophage scavenger receptor types I and II protein (271-446AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrophage scavenger receptor types I and II (Macrophage acetylated LDL receptor I and II) (Scavenger receptor class A member 1) (CD antigen CD204), MSR1, SCARA1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MSR1