



NCOA5 Antibody, HRP conjugated

Product Code	CSB-PA015544LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCD5
Immunogen	Recombinant Human Nuclear receptor coactivator 5 protein (112-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor coactivator 5 (NCoA-5) (Coactivator independent of AF-2) (CIA), NCOA5, KIAA1637
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NCOA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.